



## Gas Sorption Analyzer

**3FLEX 3500 from Micromeritics Instrument Corporation, USA** is recognized as one of the most advanced instrument in the field for material surface characterization. It is a fully automated 3-port (with 2-Micropore ports and 1-Mesopore port) gas adsorption analyzer for providing high resolution adsorption, desorption, isotherms to advance and validate fundamental understanding into development of methodologies and process improvements.

### Contact Us

Dr, A. Samanta

Department of Chemical Engineering

Phone: 326 2235087



Email: [asamanta@iitdhanbad.ac.in](mailto:asamanta@iitdhanbad.ac.in)

Web:

[https://www.iitism.ac.in/index.php/Research\\_consultancy/central\\_research\\_facility](https://www.iitism.ac.in/index.php/Research_consultancy/central_research_facility)

### Major Features:


- Minimum Measurable Surface Area Range: 0.01 m<sup>2</sup>/g
- Pore Diameter Range: 3.5 - 5000 Å
- Micropore Volume: Detectable within 0.0001 cc/g
- Analysis Range:  $1.3 \times 10^{-9}$  to 1.0 P/P<sub>0</sub>
- Ultimate Vacuum:  $3.75 \times 10^{-10}$  mm Hg



**CENTRAL RESEARCH  
FACILITY IIT(ISM)  
DHANBAD**  
Jharkhand 826004

**CENTRAL  
RESEARCH  
FACILITY  
IIT(ISM)  
DHANBAD**

*Gas Sorption Analyzer*



## Bench top Instrument is capable of carrying out following analyses:

- Single-point BET surface area
- Multi-point BET surface area
- Adsorption isotherm
- Desorption isotherms
- Langmuir surface area
- BJH mesopore volume
- BJH mesopore area
- Total pore volume
- DeBoer t-plot & others

- MP method
- Micropore volume
- Micropore area
- Micropore size distribution
- Heat of adsorption


*One instrument with three configurable analysis ports for high throughput mesopore and micropore analyses*

### Vapor Isotherms

The 3Flex includes an extensive library of fluid properties of fixed gases and commonly used vapors. Isotherm data are easily collected using hydrocarbons as the adsorptive.

### Data Reduction and Software Control

User-selectable data ranges through the graphic interface allow direct modeling for BET, BJH, t-Plot, Langmuir, DFT interpretation.



MicroActive interactive features reduce the often-difficult trial-and-error procedure for fitting the BET surface area plot by providing direct feedback to the user

