Serial No. Date:

MULTI ANGEL LIGHT SCATTERING SPECTROPHOTOMETER (MALSS) WITH GPC CENTRAL RESEARCH FACILITY INDIAN INSTITUTE OF TECHNOLOGY (ISM), DHANBAD – 826004

[Please read the information	given overleaf before filli	ng up this form and put a $()$ in appropriate box]		
I wish to get (in words) number of samples be examined/analyzed. The nature of				
sample is Water dissolve (),	, THF dissolve ().			
Please allot me slots for	samples, the estimated cha	arge for the analysis is ₹		
	Name of Supervisor			
		Cell No.:		
Signature of HOD/HOC/Guide.	/PI/Profin-charge	Signature of the user		
Dhanbad Cash counter (Receipt N	oosited through [†] DD (Numb Io) on OR	per		
† Please provide the original DD/	CC of cash receipt along w	ith this form.		
Please allot time and complete	the analysis.	Signature of the User /Faculty /Supervisor /PI		
		Signature of the Laboratory In-Charge		
The above work has been done	e satisfactorily on	(date) and generated data has been delivered		
to me.				
Signature of the Operator		Signature of the user		

INFORMATION FOR USERS

The charges for the MALSS as follows:

	MALSS analysis for Water solvent (₹)	MALSS analysis for THF solvent (₹)
For users of IIT (ISM) (per sample)	750	1000
For outside R&D* and Academics (per sample)	1000	1500
For Industry* (per sample)	2000	3000

[No GST is required for user of IIT(ISM)]

Booking Rules and Sample preparation for GPC analysis

- 1) All payment must be made prior to booking of the slot and true copy the payment slip [for deposit in IIT(ISM) cash counter in the head of CRF-MALSS] or original DD [must be drawn in favour of Registrar, IIT(ISM)] must be provided with booking form.
- 2) All forms must be forwarded through the concerned HOD, HOC, PI, Guide or Prof.-in-charge etc. and to be submitted in the MALSS Laboratory.
- 3) Polymers soluble in water/THF can be analysed.
- 4) Polymer should be completely soluble and no particle should be invisible.
 After dissolving, the polymer solution should be filtered through syringe filter (preferably 0.20 μm or less)
- 5) The solution should be as diluted as possible.
- 6) Exact concentration (upto 4 decimal) in mg/mL is required for analysis.
- 7) Polymer soluble in a particular solvent should be HPLC grade solvent e.g. if a polymer is soluble in water then HPLC grade water has to be used for dissolving. Similarly, polymer soluble in THF, HPLC grade THF has to be used.
- 8) Volume of the polymer solution should be minimum of 10 ml.

^{*}The charge are excluding GST and it may be calculated as per govt. Rule.